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## Instructions for use

### Mode d'emploi

DE

#### Intensiv Ortho-StripsTray

##### Produktbeschreibung

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- Bur-shank silicone holders orange, REF 056orange/10
- Bur-shank silicone holders yellow, REF 056yellow/10
- Bur-shank silicone holders white, REF 056white/10
- Bur-shank silicone holders red, REF 056red/10
- Bur-shank silicone holders brown, REF 056brown/10
- Bur-shank silicone holders grey, REF 056grey/10
- Bur-shank silicone holders dark blue, REF 056darkblue/10
- Bur-shank silicone holders green, REF 056green/10
- Bur-shank silicone holders light blue, REF 056lightblue/10
- Bur-shank silicone holders assorted, REF 056ass/8

EN

#### Intensiv Ortho-StripsTray

##### Product description

- Tray in stainless steel, consisting in a base and lid connected by a hinge.
- The Tray can be closed by overlapping the cover on the base.
- Once it is closed, the instruments inside are stable in their position and will not move.
- The Tray has an internal folding support with six holes.
- The holes are suitable to support sterilisable coloured silicone holders.
- The silicone holders enable to contain the instruments in a compact and steady arrangement, suitable for the safe and orderly withdrawal and relocation of the instruments.
- The colour of the silicone holders can be chosen to match the colour code of the instrument diamond grit-size.
- Dimensions:
  - Length 87mm
  - Depth 51mm
  - Height 14mm
- The Tray has slots for the disinfection liquid to flow through it.
- The Tray is designed to hold, through its silicone holders, six Intensiv Ortho-Strips, for mechanical Orthodontic stripping.

##### Intended use

- Arrangement of Intensiv Ortho-Strips during Orthodontic stripping procedures.
- Storage of Intensiv Ortho-Strips.

##### Practical tips

- Application by dental professional only.
- To open the Tray lift the cover upwards, slightly pressing the top wing and holding firmly the base.
- Move the internal folding support upwards, in order to have it parallel to the working surface. Two little wings on left and right side of the support help to tilt it in that direction.
- Place the silicone holders inside the six holes before inserting the instruments, as shown in the gures.
- Make sure that the folding internal support is moved back downwards before the Tray is closed.
- Close the Tray by pressing the cover downwards until you hear a "click".

##### Maintenance and Sterilisation

- Intensiv Ortho-StripsTray is delivered non-sterile. It must be disinfected and sterilized prior to first application to patient.
- It must be disinfected, cleaned and sterilized immediately after each use.
- If disinfection by immersion is carried out, it is appropriate to only use cleaning/disinfection solutions that provide protection against corrosion, carefully following the amounts indicated by the manufacturer.
- Sterilisation must take place according to established procedures. Use a single-pulsed or fractionated vacuum autoclave and sub-vacuum drying, according to EN 13060. Minimum exposure time at the sterilization temperature:
  - 6 minutes at 132 °C and 20 minutes drying time (US method)
  - 18 minutes at 134 °C and 20 minutes drying time (EU method)

In the unlikely case of an incident caused by the Tray please refer to the manufacturer and competent authorities.

Intensiv SA declines any liability for any use of the Tray other than that specified in these instructions for use.

Intensiv Ortho-StripsTray, REF OST400

A wide variety of different colour bur-shank silicone holders is available:

- Bur-shank silicone holders orange, REF 056orange/10
- Bur-shank silicone holders yellow, REF 056yellow/10
- Bur-shank silicone holders white, REF 056white/10
- Bur-shank silicone holders red, REF 056red/10
- Bur-shank silicone holders brown, REF 056brown/10
- Bur-shank silicone holders grey, REF 056grey/10
- Bur-shank silicone holders dark blue, REF 056darkblue/10
- Bur-shank silicone holders green, REF 056green/10
- Bur-shank silicone holders light blue, REF 056lightblue/10
- Bur-shank silicone holders assorted, REF 056ass/8

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Bur-shank silicone holders orange, REF 056orange/10  
Bur-shank silicone holders yellow, REF 056yellow/10  
Bur-shank silicone holders white, REF 056white/10  
Bur-shank silicone holders red, REF 056red/10  
Bur-shank silicone holders brown, REF 056brown/10  
Bur-shank silicone holders grey, REF 056grey/10  
Bur-shank silicone holders dark blue, REF 056darkblue/10  
Bur-shank silicone holders green, REF 056green/10  
Bur-shank silicone holders light blue, REF 056lightblue/10  
Bur-shank silicone holders assorted, REF 056ass/8

#### Bur-shank silicone holders positioning and removal